


# Hann's

## How to Measure an Accurate Dose Using an Insulin Syringe

Insulin syringes can be adjusted by either "units" or "mL." Please be sure to closely follow the conversion guidelines below to ensure proper dosage.

### MEASUREMENT CONVERSIONS

- 5 units = 0.05 mL
- 10 units = 0.1 mL
- 15 units = 0.15 mL
- 20 units = 0.2 mL
- 25 units = 0.25 mL
- 30 units = 0.3 mL



If you have any questions or need assistance, please give us a call. Our highly trained staff are available Monday through Friday 8 am to 5:30 pm and Saturday 9:00 to Noon EST.